



VOCSN Ventilation Modes

The VOCSN critical care ventilator provides a comprehensive set of ventilation modes to meet patient needs during invasive, non-invasive, and mouthpiece ventilation. VOCSN uses a standard nomenclature to describe its modes: Assist/Control (Volume or Pressure), SIMV (Volume or Pressure), Vol. Targeted (PS, PC, or SIMV), Bi-Level, and Spontaneous. The reference list below illustrates how to set your preferred ventilation mode using equivalent VOCSN modes and settings. (NOTE: This document applies to software version 4.06R and later.)

Non-Invasive Ventilation Modes

Ventilation Mode or Feature	Equivalent VOCSN Mode Settings	Additional VOCSN Controls
S/T (Spontaneous/Timed)	o Set Mode to Bi-Level o Set Breath Rate o Set Flow Trigger	o Inspiratory Time, EPAP, IPAP (measured from a baseline of 0 cmH2O), Flow Cycle, Time Cycle, and Rise Time
AVAPS TM , PRVS (Average Volume Assured Pressure Support), (Pressure Regulated Volume Support)	o Set Mode to Vol. Targeted-PS o Set Breath Rate o Set Pres. Minimum and Pres. Adj. Rate o Set Flow Trigger o Set High Pressure Alarm 5 cmH2O above desired maximum pressure	o Inspiratory Time, PEEP, Tidal Volume, Pressure Control Flow Termination, Flow Cycle, Time Cycle, and Rise Time
T (Timed)	o Set Mode to AC-Pressure o Set Breath Rate o Set Flow Trigger to Off o Set Pressure Control Flow Termination to Off	o Inspiratory Time, PEEP, Pressure Control (measured above the set PEEP), and Rise Time
Pressure Support or S (Spontaneous)	o Set Mode to SIMV-Pressure or SIMV-Volume o Set Breath Rate to 0 o Set Pressure Support	o Inspiratory Time (Apnea backup), Pressure Control (Apnea backup if SIMV-Pressure), Tidal Volume (Apnea backup if SIMV-Volume), PEEP, Flow Cycle, Time Cycle, Rise Time, Apnea Rate
CPAP (Continuous Positive Airway Pressure)	o Set Mode to SIMV-Pressure or SIMV-Volume o Set Breath Rate to 0 o Set PEEP to desired CPAP level o Set Pressure Support to 0	o Inspiratory Time (Apnea backup), Pressure Control (Apnea backup if SIMV-Pressure), Tidal Volume (Apnea backup if SIMV-Volume), Flow Cycle, Time Cycle, Rise Time, and Apnea Rate



Volume Ventilation Modes

Ventilation Mode or Feature	Equivalent VOCSN Mode Settings	Additional VOCSN Controls
AC-Volume (Assist/Control) or VC (Volume Control)	o Set Mode to AC-Volume o Set Breath Rate o Set Flow Trigger to On	o Inspiratory time, Tidal Volume, PEEP, and Sigh
SIMV-Volume or SIMV (Synchronized Intermittent Mandatory Ventilation) or	o Set Mode to SIMV-Volume o Set Breath Rate o Set Flow Trigger	o Inspiratory Time, Tidal Volume (for mandatory breaths), PEEP, Flow Cycle, Time Cycle, Rise Time, Apnea Rate, and Sigh
+/- Pressure Support	o Set Pressure Support (measured above the set PEEP)	
CV-Volume (Control Ventilation)	o Set Mode to AC-Volume o Set Breath Rate o Set Flow Trigger to Off	o Inspiratory time, Tidal Volume, PEEP, Sigh

Pressure Ventilation Modes (Including Volume-Targeted Ventilation)

Ventilation Mode or Feature	Equivalent VOCSN Mode Settings	Additional VOCSN Controls
AC-Pressure (Assist/Control) or PC (Pressure Control)	o Set Mode to AC-Pressure o Set Breath Rate o Set Flow Trigger	o Inspiratory Time, PEEP, Pressure Control (measured above the set PEEP), Pressure Control Flow Termination, Flow Cycle, and Rise Time
PRVC (Pressure Regulated Volume Control)	o Set Mode to Vol. Targeted-PC o Set Breath Rate o Set Pres. Minimum and Pres. Adj. Rate o Set Flow Trigger o Set High Pressure Alarm 5 cmH2O above desired maximum pressure	o Pres. Minimum, Pres. Adj. Rate, Inspiratory Time PEEP, Pressure Control Flow Termination, Flow Cycle, and Rise Time
SIMV-Pressure or PC-SIMV (Synchronized Intermittent Mandatory Ventilation) or +/- Pressure Support	o Set Mode to SIMV-Pressure o Set Breath Rate o Set Flow Trigger o Set Pressure Control o Set Pressure Support (measured above the set PEEP)	o Inspiratory Time, PEEP, Pressure Control Flow Termination, Flow Cycle, Time Cycle, and Rise Time
SIMV + PRVC (Synchronized Intermittent Mandatory Ventilation) + (Pressure Regulated Volume Control)	o Set Mode to Vol. Targeted-SIMV o Set Breath Rate o Set Pres. Minimum and Pres. Adj. Rate o Set Flow Trigger o Set High Pressure Alarm 5 cmH2O above desired maximum pressure	o PEEP, Tidal Volume, Pressure Support (measured above the set PEEP), Inspiratory Time, Pressure Control Flow Termination, Flow Cycle, Time Cycle, and Rise Time

VentecLife.com info@venteclife.com 1-844-MY-VOCSN



