# 510(k) Summary

Date Prepared April 6, 2017

**Sponsor** Ventec Life Systems, Inc.

19021 120th Ave. NE, Suite 101

Bothell, WA 98011 Tel: 425-686-1750

Contact Information Joseph Cipollone

VP Quality & Regulatory

19021 120th Ave. NE, Suite 101

Bothell, WA 98011 Tel: 909-213-3691

jcipollone@venteclife.com

**Device Proprietary Name:** VOCSN Unified Respiratory System

Common name Mechanical ventilator / Oxygen concentrator / Cough assist / Suction pump

Primary Classification and

**Product Code** 

Mechanical ventilator

Classification number 21 CFR 868.5895

Class II

Code: CBK - Continuous Ventilator, Facility Use

Subsequent Product Codes of Additional Integrated Functions

Mechanical ventilator

Code: NOU - Continuous Ventilator, Home Use

Oxygen concentrator

• Code: CAW - Portable oxygen generator

Cough Assist

• Code: NHJ - Noncontinuous ventilator (IPPB)

**Suction Pump** 

Code: BTA - Powered suction pump

**Heated Patient Circuits** 

Code: BZE - Heater, Breathing System, W/Wo Controller

Bacteria Filter

• Code: CAH - Filter, Bacterial, Breathing-Circuit

510(k) Summary Page 1 of 17

**Primary Predicate Device** 

Mechanical ventilator component

Information:

Respironics Trilogy Series Ventilator (K111610)

**Reference Predicate Device** 

Information:

Mechanical ventilator component

CareFusion Palmtop Ventilator PTV-8 & -10 (K070594)

CareFusion LTV 1200 Ventilator (K060647)

Oxygen concentrator component

Omni 3 (eQuinox) Oxygen System (K120785)

Cough Assist component

Respironics Cough Assist T70 (K121955)

Suction Pump component

Precision Medical Easy Go Vac PM66 (K140179)

**Heated Patient Circuits** 

Respironics (Philips) Disposable Heated Wire Circuits (K110398)

Bacteria Filter

Respironics (Philips) Bacteria Filter PN 342077 - Supplied by King

Systems under (K973797)

# Indications for Use (VOCSN Unified Respiratory System)

The VOCSN Unified Respiratory System is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive applications. The VOCSN is intended for pediatric through adult patients weighing at least 5 kg. It is intended for use in home, hospital, institutional and transport settings, including portable applications.

The integral oxygen concentrator is intended for the administration of supplemental oxygen. The integral suction pump is intended for airway fluid removal and oral/ pharyngeal hygiene. The integral cough assist option is intended for patients who are additionally unable to cough or clear secretions effectively.

#### Discussion of Differences in the Indications for Use from the Primary Predicate Device

The VOCSN and primary predicate Trilogy ventilator have substantially equivalent indications for use with regard to the ventilation therapy, except that the VOCSN is indicated for transport settings whereas the Trilogy is not intended for use as a transport ventilator. This additional use setting is equivalent to a secondary predicate that is identified.

The indications for use statement of the VOCSN encompasses additional intended uses for oxygen concentrator, cough assist, and suction pump therapies that are not available in the Trilogy ventilator. These therapies are complementary adjunct therapies to the ventilation therapy, for which secondary predicates are identified.

510(k) Summary Page 2 of 17

The Trilogy ventilator indicates an optional Oximetry Interface Kit that is not indicated in the VOCSN indications for use. Oximetry is an ancillary therapy that is not critical to the therapeutic benefit of the ventilation, oxygen concentration, cough assist and suction therapies provided by VOCSN.

# **Device Description Overview**

The VOCSN unified respiratory support system is a mechanical ventilator which combines additional conventional therapies into a single device. Additional therapies include oxygen, cough assist, and suction.

The device description will be broken down by therapy; i.e., ventilation, oxygen concentration and delivery, cough assist, and suction.

510(k) Summary Page 3 of 17

#### Ventilation

## Description

The ventilator function of the VOCSN device is a conventional positive pressure ventilator that supplies volume breaths, pressure breaths, and spontaneous breaths. These breath types are incorporated into the following traditional modes of ventilation:

- Assist Control
- SIMV (Synchronized Intermittent Mandatory Ventilation)
- Bi-Level

The VOCSN uses a conventional radial blower working in conjunction with valves and transducers under the control of a microprocessor to provide breath delivery. Oxygen is provided to the patient from the internal concentrator, or from external sources. The ventilator is connected to the patient via one of three available circuits: an active circuit, a passive circuit, or a mouthpiece circuit. An optional Secretion Trap can be placed between the patient connection and the patient circuit to collect liquids to help prevent obstructing the patient circuit. The VOCSN provides a conventional alarm and monitoring package to alert users to hazard conditions.

A central user interface provides adjustment of the controls and display of monitored data. The VOCSN is a portable device that can be operated from common AC and DC supply power sources, as well as internal batteries.

#### **Substantial Equivalence**

The intended use, performance, and technology of the ventilator system have been compared to predicate devices, the Trilogy Series Ventilator (K111610), and the Palmtop Ventilator PTV-8 & -10 (K070594). A summary table of key characteristics compared to the predicate device(s) is shown below.

Characteristic	VOCSN	Predicate, Trilogy (K111610)	Comparison
Intended Use	<ul> <li>continuous or intermittent ventilatory support</li> <li>invasive and non-invasive</li> <li>ped through adult ≥ 5 kg</li> <li>home, hospital, institutional and transport settings, including portable applications</li> </ul>	<ul> <li>continuous or intermittent ventilatory support</li> <li>invasive and non-invasive</li> <li>ped through adult ≥ 5 kg</li> <li>home, hospital, and mobile applications.</li> </ul>	Substantially Equivalent Equivalent except the Trilogy does not specify transport settings. The secondary predicate Palmtop ventilator specifies <u>transport</u> use.
FDA Product Code	NOU, CBK	NOU, CBK	Equivalent

510(k) Summary Page 4 of 17

Characteristic	VOCSN	Predicate, Trilogy (K111610)	Comparison
Modes of			Substantially Equivalent
ventilation			Facilitate
	<ul> <li>Spontaneous</li> </ul>	Spontaneous ventilation (S)	Equivalent
	Bi-Level	Spontaneous ventilation	VOCSN Bi-level mode is
		with timed back-up (S/T)	equivalent to Trilogy (S/T)
		Timed ventilation (T)	and (T) modes with
			appropriate settings of patient trigger
			patient trigger
	Assist/Control-Pressure	Pressure A/C (Palmtop	Equivalent to the secondary
		ventilator K070594)	predicate Palmtop ventilator.
	Assist/Control-Volume	A / C / D . L .	Equivalent to the secondary
	SIMV-Pressure	Volume A/C (Palmtop ventilator K070594)	predicate Palmtop ventilator.
	Sliviv-Pressure	ventuator ko70394)	predicate raintop ventuatori
	SIMV-Volume	Pressure SIMV (Palmtop	Equivalent to secondary
		ventilator K070594)	predicate Palmtop ventilator.
		Volume SIMV (Palmtop	Equivalent to secondary
		ventilator K070594)	predicate Palmtop ventilator.
		Tentinates nevees i,	
			Reference Tab 5, TPR-00049,
			TPR-00050 and TPR-00055 for test results
Significant	Breath Rate	Breath Rate	Substantially Equivalent
breath control	PEEP/ EPAP	PEEP/ EPAP	
parameters	Pressure/Pressure Control/	Pressure/Pressure Control/	
	IPAP  Inspiratory Time	IPAP     Inspiratory Time	
	<ul><li>Inspiratory Time</li><li>Sigh</li></ul>	<ul><li>Inspiratory Time</li><li>Sigh</li></ul>	
	Tidal volume	Tidal volume	
	• FIO2	• FIO2	
Core technology	Conventional radial blower	Conventional radial blower	Substantially Equivalent
	working in conjunction with	working in conjunction with	
	valves and transducers under the control of a microprocessor	valves and transducers under the control of a microprocessor	
	to provide breath delivery.	to provide breath delivery.	
Circuits	Single limb with active exh.	Single limb with active exh.	Equivalent
	valve	valve	
	Single limb with passive	Single limb with passive	
	<ul><li>exh. valve</li><li>Single limb with</li></ul>	<ul><li>exh. valve</li><li>Single limb with</li></ul>	
	mouthpiece	mouthpiece	
Circuit Interfaces	Invasive and non-invasive	Invasive and non-invasive	Equivalent
User Interface	LCD touch screen with	LCD screen with hard keys	Substantially Equivalent
	additional hard keys		The Palmtop ventilator
			(K070594) includes a touch screen
			3016611

510(k) Summary Page 5 of 17

Characteristic	VOCSN	Predicate, Trilogy (K111610)	Comparison
Nebulizer	Provides pneumatic drive to external OEM 6 L/min nebulizer	No nebulizer drive provided	Substantially Equivalent The alternate predicate
			Palmtop ventilator provides a pneumatic drive to an external OEM 6 L/min nebulizer
Power	AC, DC, and internal battery	AC, DC, and internal battery	Equivalent

The company completed Validation testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. This testing included biocompatibility and comparative waveform testing. Biocompatibility testing of breathing gases comprised evaluation of volatile organic compounds (VOC), particulate matter, carbon dioxide, carbon monoxide and ozone. Biocompatibility testing of the portion of the gas pathway in contact with humidification or aerosolized medications additionally included cytotoxicity, sensitization, intracutaneous reactivity, systemic toxicity, bacterial reverse mutation genotoxicity, and chemical characterization with risk assessment. These tests demonstrated that the VOCSN ventilator is compliant with the user and product requirements.

The product is compliant with and has been verified to the following standards:

Critical care ventilators	ISO 80601-2-12
Home care ventilators	ASTM F1246
Anesthetic and respiratory equipment Conical connectors	ISO 5356-1
Respiratory gas monitors	ISO 80601-2-55
Nebulizing system	ISO 27427
Breathing sets and connectors	ISO 5367
Medical electrical equipment	IEC 60601-1
Electromagnetic compatibility	IEC 60601-1-2
Medical devices for home use	IEC 60601-1-11
Alarm system	IEC 60601-1-8

510(k) Summary Page 6 of 17

## Oxygen Concentrator

#### Description

An optional internal oxygen concentrator can be selected to provide oxygen to the patient. Oxygen from the internal oxygen concentrator is delivered as a pulse dose via a tube directly to the patient interface.

The oxygen is separated from the nitrogen in room air using a conventional Pressure Swing Adsorption (PSA) oxygen concentration process. The system consists of a reciprocating compressor, sieve bed, and valves under the control of the microprocessor system. As part of the VOCSN system it uses the central user interface and can be operated from common AC and DC supply power sources, as well as internal batteries. The VOCSN provides an alarm and monitoring package to alert users to hazard conditions.

#### Substantial Equivalence

The intended use, performance, and technology of the oxygen concentrator system have been compared to the Omni 3 (eQuinox) Oxygen System (K120785) predicate device. A summary table of key characteristics compared to the predicate device is shown below.

Characteristic	VOCSN	Predicate, Omni 3 (K120785)
Intended Use	The integral oxygen	The administration of
	concentrator is intended for the	supplemental oxygen.
	administration of supplemental	
	oxygen	
FDA Product	CAW	CAW
Code		
Oxygen %	Nominal 90 %	Nominal 90 %
Modes	Pulse Dose	Pulse Dose
		<ul> <li>Continuous</li> </ul>
Core technology	Pressure Swing Adsorption (PSA)	Pressure Swing Adsorption
	oxygen concentration process.	(PSA) oxygen concentration
		process.
User Interface	LCD touch screen with	LCD screen with hard keys
	additional hard keys	
Power	AC, DC, and internal battery	AC, DC, and internal battery

#### Performance Testing (non-clinical)

The company completed Validation testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. This testing included biocompatibility testing, as listed in the Ventilation section. These tests demonstrated that the VOCSN oxygen concentrator is compliant with the user and product requirements.

The oxygen concentrator option is compliant with and has been tested to the following standards:

Oxygen concentrators

ISO 80601-2-69

510(k) Summary Page 7 of 17

Oxygen conserving equipment	ISO 80601-2-67
Medical electrical equipment	IEC 60601-1
Electromagnetic compatibility	IEC 60601-1-2
Medical devices for home use	IEC 60601-1-11
Alarm system	IEC 60601-1-8

510(k) Summary Page 8 of 17

# **Cough Assist**

# Description

The Cough Assist option clears secretions from the lungs by applying positive pressure (insufflation) to the airway followed by a sudden negative pressure (exsufflation). This creates a high expiratory flow, simulating a natural cough. The device attaches to standard patient interfaces using the VOCSN patient circuits.

The VOCSN uses a conventional radial blower working in conjunction with a cough assist valve under the control of a microprocessor to provide the insufflation and exsufflation cough phases. The cough assist function is integral to the ventilator and is connected to the patient via the patient circuit.

The VOCSN provides an alarm and monitoring package to alert users to hazard conditions.

## **Substantial Equivalence**

The intended use, performance, and technology of the cough assist system have been compared to the Cough Assist T70 (K121955) predicate device. A summary table of key characteristics compared to the predicate device is shown below.

Characteristic	VOCSN	Cough Assist T70 (K121955)
Intended Use	Mechanically ventilated patient unable to cough or clear secretions effectively	patient unable to cough or clear secretions effectively.
	<ul> <li>pediatric through adult ≥ 5 kg</li> <li>invasive and non-invasive</li> </ul>	<ul> <li>adult or pediatric patients</li> <li>used with facemask, mouthpiece, endotracheal,</li> </ul>
	<ul> <li>home, hospital, institutional</li> </ul>	<ul> <li>or tracheostomy tube.</li> <li>hospital, institutional environment, or in the</li> </ul>
	and transport settings, including portable applications	home.
FDA Product Code	NHJ	NHJ
Significant control parameters	<ul><li>Breath Sync</li><li>Exsufflation Pressure</li><li>Exsufflation Time</li></ul>	<ul> <li>Cough Trak</li> <li>Exsufflation Pressure</li> <li>Exsufflation Time</li> </ul>
	<ul> <li>Insufflation Rise Time</li> <li>Insufflation Time</li> </ul>	Inhale Flow (controls rise time)     Insufflation Time
	Pause Time	Pause Time

510(k) Summary Page 9 of 17

Characteristic	VOCSN	Cough Assist T70 (K121955)
Core technology	Radial blower working in	Radial blower working in
	conjunction with a cough assist	conjunction with two cough
	valve under the control of a	assist valves under the control
	microprocessor	of a microprocessor
User Interface	LCD touch screen with	LCD screen with hard keys
	additional hard keys	
Power	AC, DC, and internal battery	AC, DC, and internal battery

The company completed validation testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. This testing included biocompatibility testing, as listed in the Ventilation section. These tests demonstrated that the VOCSN is compliant with the user and product requirements.

The cough assist option is compliant with and has been verified to the following standards:

Medical electrical equipment IEC 60601-1

Medical devices for home use IEC 60601-1-11

Alarm system IEC 60601-1-8

510(k) Summary Page 10 of 17

#### Suction

#### Description

The VOCSN includes an optional integrated suction function and detachable VLS suction canister. If enabled, Suction therapy can be initiated at any time during Ventilation or Cough therapy, to help remove secretions from the patient airway or patient circuit.

When enabled, the Cough+Suction feature activates both therapies simultaneously to help remove secretions from the patient airway and/or patient circuit. Suction will begin at the start of the first Cough therapy insufflation, and run throughout the configured number of Cough Cycles plus an additional 30 seconds.

The VOCSN suction function uses a conventional reciprocating piston pump working in conjunction with a selector valve under the control of a microprocessor to provide the negative pressure. Suction is routed to the canister, and a suction tube connects the suction catheter to the canister.

The VOCSN provides an alarm and monitoring package to alert users to hazard conditions.

## **Substantial Equivalence**

The intended use, performance, and technology of the suction function have been compared to the Precision Medical Easy Go Vac PM66 (K140179) predicate device. A summary table of key characteristics compared to the predicate device is shown below.

Characteristic	VOCSN	Easy Go Vac PM66 (K140179)
Intended Use	<ul> <li>intended for airway fluid removal and oral/pharyngeal hygiene.</li> <li>home, hospital, institutional and transport settings, including portable applications</li> </ul>	<ul> <li>provides a portable, AC/DC powered medical vacuum source.</li> <li>It is intended for use in the homecare / healthcare environments</li> </ul>
FDA Product Code	ВТА	ВТА
Vacuum Control	-50 to -300 mmHg	-50 to -533 mmHg
Core technology	reciprocating piston pump working in conjunction with a selector valve under the control of a microprocessor	reciprocating piston pump working in conjunction with a vacuum regulator and mechanical gage
User Interface	LCD touch screen with additional hard keys	Manual physical controls
Power	AC, DC, and internal battery	AC, DC, and internal battery

510(k) Summary Page 11 of 17

The company completed Validation testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. These tests demonstrated that the VOCSN is compliant with the user and product requirements.

The suction option is compliant with and has been tested to the following standards:

FDA Guidance Document for Powered Suction Pump 510(k)s

Electrically powered suction equipment

Medical electrical equipment

Medical devices for home use

Alarm system

IEC 60601-1-8

510(k) Summary Page 12 of 17

# **Heated Wire Patient Circuit Accessory**

## Description

The VOCSN system includes a heated patient circuit accessory to connect the ventilator to the patient connection. The circuits are single patient use pediatric & adult heated patient circuits, with a passive or active exhalation valve. Passive circuits use a fixed leak for exhalation. Active circuits use a piloted diaphragm valve for exhalation control. The VOCSN heated wire patient circuits are to be used only with the VOCSN Unified Respiratory System.

Operating principle  Connects ventilator to patient connection. Passive circuit uses a fixed leak for exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.  Classification  Product Code  Intended use  The patient circuit is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive applications. The  Connects ventilator to patient connection. Passive circuit uses a fixed leak for exhalation connection. Passive circuit uses a fixed leak for exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.  Class II  Product Code  BZE  Intended use  The patient circuit is intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive ventilation			Respironics (Philips) Disposable
principle    connection. Passive circuit uses a fixed leak for exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.    Classification   Class II   Class II   Product Code   BZE   Intended use   The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive   connection. Passive circuit uses fixed leak for exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.    Class II   Class II   Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive			Heated Wire Circuits (K110398)
uses a fixed leak for exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.  Classification Class II Class II  Product Code BZE BZE  Intended use The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive			•
exhalation port. Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.  Classification  Class II  Product Code  Intended use  The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  Active circuit uses a piloted diaphragm valve for exhalation control. Includes heated wire for exhalation control. Includes heated wire for use with a humidifier.  Class II  Product Code  BZE  BZE  Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive	principle		
uses a piloted diaphragm valve for exhalation control. Includes heated wire for use with a humidifier.  Classification Class II Class II  Product Code BZE BZE  Intended use The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive			•
valve for exhalation control. Includes heated wire for use with a humidifier.  Classification Class II  Product Code BZE Intended use The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  control. Includes heated wire for use with a humidifier.  Class II  Product Code BZE Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive		-	
Includes heated wire for use with a humidifier.  Classification Class II Class II  Product Code BZE BZE  Intended use The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive used invasive and non-invasive in termittent ventilation. It may be used for both invasive and non-invasive			. •
with a humidifier.  Class II  Product Code  BZE  Intended use  The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  With a humidifier.  Class II  Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive		valve for exhalation control.	control. Includes heated wire for
Classification  Product Code  BZE  Intended use  The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  Class II  BZE  Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive			use with a humidifier.
Product Code  Intended use  The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  The patient circuit is intended to and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive	Classification		Class II
Intended use  The patient circuit is intended to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  Intended to provide warmed and/or humidified breathing gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive			
to be used only with the VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non- individuals who reads in invasive and non- invasive and non- invasive and non- invasive and non- invasive			
VOCSN Unified Respiratory System which is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive  gases before they enter the patient's airway. The disposable heated wire circuit is indicated for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive	interiaca asc	•	•
System which is intended to provide continuous or heated wire circuit is indicated for use by a single adult or support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive		•	
provide continuous or intermittent ventilatory for use by a single adult or support for the care of individuals who require mechanical ventilation. It may be used in invasive and non-invasive		• • •	•
intermittent ventilatory support for the care of individuals who require mechanical ventilation. It may be used in invasive and non- intermittent ventilatory for use by a single adult or pediatric patient in the home, hospital, and or institutional setting. It may be used for both invasive and non-invasive		= -	
support for the care of pediatric patient in the home, individuals who require hospital, and or institutional mechanical ventilation. It may be used for both be used in invasive and non-invasive		•	
individuals who require hospital, and or institutional mechanical ventilation. It may be used for both be used in invasive and non-invasive		,	
mechanical ventilation. It may setting. It may be used for both be used in invasive and non-invasive			
be used in invasive and non-invasive		•	•
			,
invasive applications. The ventuation			
VOCSN is intended for			ventuation
pediatric			
through adult patients		•	
weighing at least 5 kg. It is		= -	
intended for use in home,			
hospital, institutional and		*	
transport settings, including			
portable applications.			
Humidifier Fisher Paykel HC550 & HC500	Humidifier		HC500
compatibility MR850 humidifiers		•	116566
Reuseable Single patient use Single patient use			Single patient use
	-	= -	Applied voltage through heating
heating wires. wires.			
Sterility Non-Sterile Non-Sterile	Sterility		Non-Sterile
Materials Compliant with ISO 10993-1 Compliant with ISO 10993-1			
Biological Evaluation Biological Evaluation		•	•
Breathing gases	Breathing gases	=	=

510(k) Summary Page 13 of 17

Feature	VOCSN	Respironics (Philips) Disposable Heated Wire Circuits (K110398)
Tube diameter	Adult: 22 MM Pediatric 15 mm	Adult: 22 MM Pediatric 15 mm
Tube length	1.85 M	1.8 m
Tube	22 mm conical Compliant	22 mm conical Compliant with
connectors	with ISO-5356-1 Conical connectors	ISO-5356-1 Conical connectors
Heating wire	Resistive, encased	Resistive, encased
Wire resistance	16	30
(ohms)		
Rated power	60 watts	Not stated
Patient Leakage	Compliant with the	Claims compliance with IEC
current	requirements of IEC 60601-1	60601-1
Resistance to airflow	Complaint with ISO 8185	Complaint with ISO 5367
	Pediatric: approx. 326 ml	Not stated, claims compliance
Tube volume	· ·	•
	Adult: approx. 703 ml	with ISO 5367
	Compliant with ISO 5367	
Exhalation valve	Passive Option – fixed orifice	Passive Option – fixed orifice type
type	type	A .: O .: SII . I
		Active Option – Piloted
	Active Option – Piloted	diaphragm type
	diaphragm type	
Enthalpy	Per ISO 8185	Not stated
Condensate	The VOCSN circuit was	Not stated.
performance	demonstrated to control	
	condensate over a wide range	
	of patient ventilation	
	conditions.	
ISO-5367	Complies	Complies
Breathing tubes		
ISO-8185	Complies as applied to	Complies as applied to breathing
Respiratory	breathing circuits	circuits
humidification		
systems		
ISO-60601-1	Complies with applicable	Complies with applicable parts
Medical	parts as specified in ISO 8185	
Electrical		
Equipment		

The company completed Validation testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. This testing included biocompatibility testing, as listed in the Ventilation section. These tests demonstrated that the VOCSN is compliant with the user and product requirements.

510(k) Summary Page 14 of 17

The heated wire patient circuit accessory is compliant with and has been tested to the following standards:

Biological Evaluation	ISO 10993-1
Respiratory tract humidifiers for medical use	ISO 8185
Conical connectors	ISO 5356-1
Breathing sets and connectors	ISO 5367

510(k) Summary Page 15 of 17

# **Bacteria Filter Accessory**

## Description

The VOCSN system includes a bacteria filter accessory connected between the VOCSN and the patient circuit. The filter is intended to reduce bacterial/viral transmissions between the patient and equipment. The VOCSN bacteria filter is to be used only with the VOCSN Unified Respiratory System.

Feature	VOCSN	Respironics (Philips) Bacteria Filter PN 342077 - Supplied by King Systems under (K973797)
Classification	Class II	Class II
Product Code	CAH	CAH
Intended use	To be used only as part of the VOCSN integrated respiratory care system to reduce bacterial/viral transmissions between the patient and equipment.	It is designed to reduce bacterial/viral transmissions between the patient and equipment.
Operating principle	Electrostatic	Electrostatic
Filtration efficiency	99.99% BFE & VFE	99.99% BFE & VFE
Approximate Volume (ml)	31	74
Connectors	22 mm conical, compliant with ISO 5356-1 Conical Connectors	22 mm conical
Resistance @ 30 lpm	0.4 cmH2O	0.7 cmH2O

### Performance Testing (non-clinical)

The company completed testing to validate performance against the User Requirements Document under simulated use conditions, and test reports were generated summarizing the test results. Verification testing was completed to verify system performance against the Product Requirements Specification, and test reports were generated summarizing the test results. This testing included biocompatibility testing, as listed in the Ventilation section. These tests demonstrated that the VOCSN is compliant with the user and product requirements.

The bacteria filter accessory is compliant with and has been tested to the following standards:

Conical connectors ISO 5356-1

510(k) Summary Page 16 of 17

## Conclusion

The indications for use, performance characteristics, and core technology of the VOCSN are substantially equivalent to the predicate devices. Extensive validation and verification testing has demonstrated that the device is compliant with the product requirements and relevant regulatory standards. The information provided supports the claim that the device is substantially equivalent to predicate devices.

510(k) Summary Page 17 of 17