

# HOME VENTILATORS

## Featured Products

\*HCPCS codes subject to change



1




2



3



4



6

### 1 VENTEC LIFE SYSTEMS

#### VOCSN

VOCSN integrates five separate medical devices, including a critical care ventilator, 6L/min equivalent oxygen concentrator, touch-button cough assist, hospital grade suction and a high-performance nebulizer, into one unified respiratory system.

Lightweight for easy transport, VOCSN is controlled with touchscreen functionality to switch between therapies. Coded HCPCS E0467. Visit [venteclife.com](http://venteclife.com).

Circle 203 on card.

### 2 BREAS

#### Vivo 65

The Vivo 65, an advanced life-support ventilator, combines clinical excellence, connectivity and monitoring possibilities for both pediatric (>11 lbs / 5 kg) and adult patients in hospitals, post-acute care institutions and the home. The Vivo 65's combined integrated monitoring of SpO<sub>2</sub>, EtCO<sub>2</sub> and FiO<sub>2</sub> reduces the need for additional devices. The user interface makes the device easy to use and easy to teach. Coded HCPCS E0465 and E0466.

Visit [breas.com](http://breas.com).  
Circle 204 on card.

### 3 RESMED

#### Astral

ResMed's Astral life support ventilator offers mobility and ease of use for patients suffering from neuromuscular disease, COPD and other breathing-related diseases. It weighs just 7.1 pounds with an eight-hour internal battery and 24-hour battery life with external battery connections. The Astral also has AutoEPAP, easy mouthpiece interchange from a single limb circuit, and cloud connectivity for remote patient monitoring via ResMed's AirView platform, enabling physicians and HMEs to access patient data in the cloud to better manage and improve outcomes. Coded HCPCS E0465 and E0466. Visit [resmed.com](http://resmed.com).

Circle 205 on card.

### 4 PHILIPS RESPIRONICS

#### Trilogy

Philips connected Trilogy is a portable ventilation solution that provides invasive and non-invasive ventilator support for a wide range of adult and pediatric patients in the home and at alternative care sites. Connected to Philips' cloud-based Care Orchestrator patient management application, connected Trilogy turns medical-grade data into actionable information. This solution enables care teams to coordinate and monitor patients remotely, allowing for fast and informed clinical decisions, including early intervention, which can help avoid unnecessary readmissions and lower cost of care. Coded HCPCS E0465 and E0466. Visit [usa.philips.com](http://usa.philips.com).

Circle 206 on card.

### 5 BREATHE TECHNOLOGIES

#### Life2000 Ventilation System

The Life2000 Ventilation System is a modular ventilation system that provides ventilatory support while facilitating ambulation and activity for up to 24 hours of use. Life2000 was designed to allow patients to participate more fully in life and perform activities of daily living (ADLs) inside and outside the home. The one-pound, wearable Life2000 Ventilator may alleviate COPD-related symptom burden and provide patients with the opportunity to perform more activity. The low-profile nasal pillow interface is similar in style and tubing size to what the patient is accustomed to using when compared to their oxygen nasal cannula. Coded HCPCS E0466 and E0465. Visit [breathetechnologies.com](http://breathetechnologies.com).

Circle 207 on card.

### 6 FISHER & PAYKEL HEALTHCARE

#### myAIRVO

myAIRVO humidified high flow combines humidification and respiratory support for spontaneously breathing patients that could benefit from high-flow therapy and airway hydration. By providing optimal humidification at 37°C/44mg/L, and respiratory support from 10-60 L/min of air flow, myAIRVO aids in secretion mobilization while meeting patients' peak inspiratory demand to slow their respiratory rate and improve alveolar ventilation. myAIRVO delivers comfortable, warm air via nasal cannula or tracheostomy interface, and is a compact unit that can be used at bedside. Supplemental oxygen can be administered with an oxygen concentrator, when necessary. Visit [fphcare.com/us](http://fphcare.com/us).

Circle 208 on card.